



## Disposable Forehead Thermometer Strips

- Liquid Crystal Technology Thermometer with Medical Grade Adhesive
- All Paper Thermometer - Mylar Free
- No Mercury, No Glass, No Batteries and Latex-Free
- 94° F - 104° F Range - Accurate to Within 1 Degree
- Designed for Single Use to Minimize Contamination Risk
- Continuous Monitoring Overnight
- Private Labeling Available for Retail Environments or Promotional Use
- Alternate Packaging and Display Options Available



**Now Taking Pre-Orders**  
(Please allow 3-4 weeks lead time for shipment)

**Call now for pricing!**  
**813.962.8890**



# EASY-TO READ, REUSABLE FOREHEAD THERMOMETERS

## INSTRUCTIONS FOR F1001 REUSABLE MYLAR THERMOMETER STRIPS

Thermometers should be used indoor and at room temperature (approximately 68° - 86°F). Avoid direct sunlight or strong heat when taking temperature. Do not take temperature if individual has been eating, drinking or exercising within the last 30 minutes due to the possibility of temperature fluctuation.

- 1. Hold thermometer at both ends, pressing firmly on dry forehead for approximately 15 seconds.**
- 2. Read thermometer while on forehead:**
  - **Green: correct temperature**
  - **If you do not see green, blue or tan colors will appear:**
    - Blue: add 1°F**
    - Tan: subtract 1°F**

Always consult your physician if a fever exists.  
The reusable thermometer may be cleaned with a soft cloth.

**TO BE ADMINISTERED BY ADULTS ONLY.  
KEEP OUT OF REACH OF CHILDREN.**

MANUFACTURED IN THE U.S.

